



Disposable Safety-Cap Needle(Syringe)



Saferwith

Provides safety not only medical but also other fields.

Safety Engineered Injection Devices (so-called 'Safety Syringes')

Safety engineered injection devices["Auto Disable" (AD) and "Reuse Prevention" (RUP) syringes] specifically designed for immunization and for therapeutic injections have been developed and identified as one of the strategies to avoid syringes and needles reuse.


These syringes can also have **Sharp Injury Protection mechanisms (RUP/SIP syringes)** for the prevention of needle-stick injury among health-care workers.

(WHO, Injection Safety Policy and Global Campaign, 2015)

Disposable Safety-Cap Needle(Syringe)

Only enhanced needle cap makes it possible to work as safety syringe.

→ This **"SMART" safety cap** will make all injections safer.




Operation Procedures of Disposable Safety-Cap Needle(Syringe)



Disposable Safety-Cap Needle(Syringe) : 2MW Series

Advantages


- 1. Simplicity** : Operated in the same manner as regular syringe
- 2. Compatibility** : Compatible with almost all syringe with needle
- 3. Dual Safety System** : Provides both RUP and SIP function



4. Safer Design for SIP

Prevents NSIs(Needle-Stick Injuries), especially during recapping

- The grip position located in the middle of safety-cap prevents HCWs hands from accessing to the needle point during the whole injection process



Comparison of Recapping Process



Regular Syringe

VS



Safety-Cap Syringe


5. Multi-Usage Available (Product Simplification)

- Price and investment competitiveness



6. Presence of an Indicator

- Eliminate potential risks for user convenience



Disposable Safety-Cap Needle(Syringe) : 2MW Series

7. One Handed Opening/Recapping Devices - User convenience & Safer SIP function



Products

1. Safety-Cap Needle



2. Safety-Cap Syringe



Spec

Model	Gauge														Length		Syringe(ml)	
	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	Equipped		Available
2MW1 Series	2MWL1			√	√	√	√	√	√	√	√					60mm	35.8mm 32mm	1 / 2 / 2.5 3 / 5 / 10 20 Lure Slip
	2MWM1	√	√	√	√	√	√	√	√	√	√	√	√	√	√	38mm	25mm 20mm	
	2MWS1		√	√	√	√	√	√	√	√	√	√	√	√	√	25mm	17.4mm 16mm 13mm 8mm	
2MW3 Series	2MWS3	√	√	√	√	√	√	√	√	√	√	√	√	√	√	25mm	25mm 20mm 16mm 13mm 8mm	1 / 2 / 2.5 3 / 5 / 10 20 Lure Slip Lure Lock
	2MWI3 (Insulin)													√	√	8mm	8mm	



Room 305, BI Center, Kyonggi University, 154-42, Gwanggyosan-ro, Yeongtong-gu, Suwon-si, Gyeonggi-do, South Korea
 TEL : +82-70-8839-2850 FAX : +82-70-8839-2851 E-mail : ktpark@saferwith.com Homepage : www.saferwith.com

Homepage QR Code

